



BioflowTM *Slim*

Touch-free Liquid/Gel Dispenser





Touch-free Liquid/Gel Dispenser

Bioflow™ is a state of art, touch-free and wall mounted dispenser to dispense handrub / handwash in medical and industrial settings.

Bioflow™ automatically dispenses both liquids and gels at a prefixed dose. This ensures adequate disinfection of hands without contaminating the environment.

Bioflow™ Product Description

Features	Benefits
Touch-free	Prevents cross contamination.
After sales service	Peace of mind.
Stainless Steel	Corrosion resistance, good strength, long term value, fire & heat resistance.
Fixed dose dispensing	Adequate disinfection. Reduced wastage of handrub/ handwash.
AC adaptor provided	Battery optional.
Compatible with liquids and gels	Versatile.

Product Specifications	
Product Name	: Bioflow™ Slim
CAT No.	: 844BS000003
Unit Size	: 173 L x 90 W x 290 H (mm)
Weight	: 2 kg
Material	: SS 304-Assembly
Bottle Volume	: 500 ml
Installation	: Wall mount / Table top
Power Supply	: Input 100-270V AC & Output 5V, 2 AMP DC
Battery (Li-Ion)	: 3.7 V (2200 mAh)

* Only Machine to provide disposable bottle & nozzle assembly.



BIOFLOW™ is compatible with

ALCONOX®

Colourless & odourless alcoholic handrub with moisturiser

ECOMAX™

Alcoholic handrub with moisturiser

PURELLIUM™ GEL

Alcoholic handrub gel with moisturiser

STERIMAX®

Liquid handrub antiseptic with triple action

TRIOSEPT™

Colourless & odourless liquid handrub antiseptic with triple action

BIOSCRUB™

Antiseptic surgical scrub

HITMAX™

Liquid microbicidal handwash soap

A Division of Tulip Diagnostics (P) Ltd.



Applying Science In Disinfection

Gitanjali, Tulip Block, Dr. Antonio do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403 202, INDIA.
Tel.: +91 832 2458546-50 Fax : +91 832 2458544 E-mail : sales@tulipgroup.com Website : www.tulipgroup.com